

## Creatinine (Jaffe modified method)

In-vitro Diagnostic reagent/kit for quantitative determination of Creatinine in serum/plasma and urine sample on Photometric System.

### Reagent

Reagent 1: Alkaline Solution  
 Reagent 2: Picric Acid Solution  
 Standard: Creatinine (Conc. 2 mg/dL)

### Principle

Creatinine forms a colored complex with picrate in alkaline medium. The rate of formation of the complex is measured photometrically.

### Summary

Creatinine is filtered by kidneys as waste product. Thus the concentration of creatinine in blood/serum of a normal individual is fairly constant. Therefore, increased blood/serum creatinine values always indicate decreased excretion meaning impaired kidney function. The creatinine concentration enables a quite good estimation for the detection of kidney diseases and monitoring of renal function. For this purpose creatinine is measured simultaneously in serum and urine collected over a defined time period.

### Storage Instructions and Reagent Stability

The reagents and standard are stable till the date of expiry, if stored at 2°C-30°C, protected from light and contamination is avoided. Do not freeze the reagents.

### Composition

Reagent – Picric acid - 9 mmol/L, Sodium hydroxide - 0.4 mol/L.  
 Standard – Creatinine 2 mg/100mL

### Waste Management

Please refer to local regulatory requirements.

### Reagent Preparation

Mix 1 part of reagent 1 (R1) + 1 part reagent 2 (R2) = Working reagent. The working reagent is stable for 5 days at 2°C–25°C when stored in dark bottle.

The working reagent must be protected from light.

### Materials required but not provided

NaCl solution 9 g/L  
 General laboratory equipment

### Specimen

Serum, heparin plasma or EDTA plasma Stability:

In Serum/Plasma 7 months at 2° – 8°C,  
 3 months at -20°C  
 In Urine 6 Days at 2° – 8°C,  
 3 months at -20°C

Dilute Urine 1 + 49 with distilled water  
 Discard contaminated specimens.

### Assay Procedure

Wavelength 505 nm (490-510 nm)  
 Optical path 1 cm  
 Temperature 37°C

Sample/Standard	100 µL
Working reagent	1000 µL
Mix, incubate for 20 sec and read absorbance (A1) & exactly after another 60 sec (A2)	

### Calculation

Calculation of the concentration “C” of Creatinine in serum or plasma.

$$C = 2.0 \times \frac{\Delta A \text{ Sample}}{\Delta A \text{ Standard}} \text{ [mg/dL]}$$

Calculation of the concentration “C” of Creatinine in urine.

$$C = 2.0 \times 50 \times \frac{\Delta A \text{ Sample}}{\Delta A \text{ Standard}} \text{ [mg/dL]}$$

$$\text{Creatinine clearance} = \frac{(\text{mg creatinine/dL urine}) \times (\text{mL urine/24 hr.})}{(\text{mg creatinine/dL serum}) \times 1440} \text{ [mL/min]}$$

### Quality Control

For internal quality normal and abnormal controls should be assayed with each batch of samples.

Each laboratory should establish corrective action in case of deviations in control recovery.

### Warnings and Precautions

- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Wear suitable gloves and eye/face protection.
- Always use safety pipettes to pull the reagents into a pipette.
- Reagents may contain some non-reactive and preservative components. It is suggested to handle carefully, avoid direct contact with skin and do not swallow.
- The reagents contain sodium hydroxide. Do not swallow. Avoid contact with skin and mucous membranes.
- For professional use only!

### Performance Characteristics Measuring Range

The test has been developed to determine Creatinine activity from 0.20 mg/dL to 25 mg/dL. If such value is exceeded the sample should be diluted 1+9 with NaCl solution (9 g/L) and the result is multiplied by 10.

### Linearity/Limit of Maximum Detection

The maximum limit of detection is 25 mg/dL.

### Sensitivity/Limit of Detection

The lower limit of detection is 0.20 mg/dL.

### Interferences

No Interferences was observed by Ascorbic acid up to 30 mg/dL, Bilirubin up to 4 mg/dL and Triglycerides up to 2000 mg/dL.

### Precision

Intra assay n = 20	Mean [mg/dL]	SD [mg/dL]	CV [%]
Sample 1	0.98	0.03	3.01
Sample 2	2.02	0.05	2.6
Sample 3	5.53	0.09	1.68

Inter assay n = 20	Mean [mg/dL]	SD [mg/dL]	CV [%]
Sample 1	1.04	0.04	3.67
Sample 2	2.20	0.04	1.86
Sample 3	6.30	0.06	0.97

### Method Comparison

A comparison of Precision Biomed Creatinine (y) with a commercially available test (x) using 20 samples gave following results:  
 $y = 0.960x + 0.031$ ;  $r^2 = 0.991$ .

### Reference Range

Unit	mg/dL	μmol/L
Serum		
Men	0.7 - 1.4	53 - 97
Women	0.6 - 1.2	44 - 80
Urine		
Men	1 - 2 g / 24 hrs.	(8.84 - 17.7) mmol/24 hrs.
Women	0.8 - 1.8 g / 24 hrs.	(7.07 - 15.9) mmol/24 hrs.

### for Creatinine Clearance

Unit	mL/min	mL/sec
Men	98 - 156	1.63 - 2.60
Women	95 - 160	1.58 - 2.67

Note: It is recommended that each laboratory should establish its own reference range based on the patient population.

### Quick Reference

Parameter	Creatinine
Mode	Fixed time
Reaction slope	Increasing
Wavelength	505 nm
Path length	10 mm
Temperature	37°C
Standard conc.	2 mg/dL
Working Reagent	1 Part Reagent 1 (500 μL) 1 Part Reagent 2 (500 μL)
Working Reagent volume	1000 μL
Sample volume	100 μL
Delay	20 Sec
Rate	60 sec
No. of readings	1 Nos.
Normal range	Men 0.7-1.4mg/dL Women 0.6 - 1.2 mg/dL Urine - 1 - 2 g/24 hours
Linearity	25 mg/dL
Sensitivity	0.20 mg/dL


### Pack Size

Cat No.	Configuration	Pack
CRE00100	Reagent 1- 1 x 50mL Reagent 2 - 1 x 50mL Standard - 1 x 2mL	100mL
CRE01000	Reagent 1 - 1 x 500mL Reagent 2 - 1 x 500mL Standard - 1 x 4mL	1000mL

### Literature

- Sarre, H. (1959) Nierenkrankheiten. Georg ThiemeVerlag, Stuttgart
- Schirmeiste.J.et.al.(1964).Dtsch.Med.Wschr.89:1640
- Mazzachi BC, Peake MJ, Ehrhardt V. Reference Range and Method Comparison Studies for Enzymatic and Jaffé Creatine Assays in Plasma and Serum and Early Morning Urine. Clin. Lab. 2000; 46: 53-55.
- Swanson AF, Swartzentruber M, Nolen PA et al. Multicenter Evaluation of the Boehringer Mannheim Compensated, Rate-Blanked Creatinine/Jaffe Application on BM/Hitachi Systems. Advances in Clinical Diagnostics. 1993. Boehringer Mannheim Corporation.
- Guder WG, Zawta B. Recommendations of the Working group on Preanalytical Quality of the German Society for Clinical Chemistry and the German Society for Laboratory Medicine: The Quality of Diagnostic Samples. 1st ed Darmstadt: GIT Verlag 2001; p. 24-5,50-1.
- Levey AS, Coresh J, Greene T, Marsh J et al: Expressing the Modification of Diet in Renal Disease Study Equation for Estimating Glomerular Filtration Rate with Standardized Serum Creatinine Values. Clin Chem 2007; 53 (4): 766-72.

Version : CRE/01



Manufactured in India by:

**Precision Biomed Pvt. Ltd.**

Plot No. - 193, "Silver Soil Industrial Park", Village - Anantpura-Chimanpura,

Teh.-Chomu, District- Jaipur - 303702 (Rajasthan) India.

Ph.: +91-9982806050 Customer Care : +91-7820806050

E-mail: info@precisionbiomed.in, Website: www.precisionbiomed.in

Mfg.Lic.No.MFG/IVD/2021/000058