


PRECISION
Be Sure ...

CE


PreciCHEM-300

Semi-auto Biochemistry Analyzer

- Ideal Companion of clinical laboratories and hospitals
- High accuracy of results
- Auto sleep mode for the lamp to prolong its life
- 48 direct test access keys for user friendly operations



- 48 direct keys to access test directly from keyboard
- Lamp saving mode indicator on keyboard
- Key to Rinse flow cell at any stage of test



High resolution graphic thermal printer for multiform reports



- Realtime absorbance curve, Realtime clock, Realtime temperature display on screen
- QC program for 3 controls per tests with westguard rules, Levey-Jennings chart
- Unique feature of Kinetic graph storage facility
- Calibration graph storage for all tests
- Facility to create profiles of various test



Software calibrated peristaltic pump for less reagent consumption



- Powerful optical assembly with autosensing of cuvette or flowcell for proper positioning
- Use of Japanese photodiode and European interference filters



- USB Port for direct attachment of external printer
- RS 232 port for data transfer
- LAN Connectivity (optional)